

Website: www.upingBio.com

SC6A3 Polyclonal Antibody

| Catalog No | YP-Ab-06215 |
|--------------------|--|
| lsotype | lgG |
| Reactivity | Human;Mouse;Rat |
| Applications | WB;ELISA |
| Gene Name | SLC6A3 DAT1 |
| Protein Name | Sodium-dependent dopamine transporter (DA transporter) (DAT) (Solute carrier family 6 member 3) |
| Immunogen | Synthesized peptide derived from part region of human protein |
| Specificity | SC6A3 Polyclonal Antibody detects endogenous levels of protein. |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution | WB 1:500-2000 ELISA 1:5000-20000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | 68kD |
| Cell Pathway | Cell membrane ; Multi-pass membrane protein . Cell projection, neuron projection . Cell projection, axon . Localizes to neurite tips in neuronal cells (By similarity). Colocalizes with SEPTIN4 at axon terminals, especially at the varicosities (By similarity) |
| Tissue Specificity | Highly expressed in substantia nigra (PubMed:7637582). Expressed in axonal varicosities in dopaminergic nerve terminals (at protein level) (PubMed:17296554). Expressed in the striatum (at protein level) (PubMed:17296554). |
| Function | function:Amine transporter. Terminates the action of dopamine by its high affinity sodium-dependent reuptake into presynaptic terminals.,miscellaneous:This protein is the target of psychomotor stimulants such as amphetamines or cocaine.,online information:Dopamine transporter entry,similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,subunit:Homooligomer; disulfide-linked. Interacts with PRKCABP and TGFB1I1., |
| Background | This gene encodes a dopamine transporter which is a member of the sodium- and chloride-dependent neurotransmitter transporter family. The 3' UTR of this gene contains a 40 bp tandem repeat, referred to as a variable number tandem repeat or VNTR, which can be present in 3 to 11 copies. Variation in the number of repeats is associated with idiopathic epilepsy, attention-deficit |



UpingBio technology Co.,Ltd

😮 Tel: 400-999-8863 📼 Emall:Upingbio.163.com

Website: www.upingBio.com

| | hyperactivity disorder, dependence on alcohol and cocaine, susceptibility to Parkinson disease and protection against nicotine dependence.[provided by RefSeq, Nov 2009], |
|---------------------------|--|
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |

Products Images