




ABHD14A Polyclonal Antibody

Catalog No	YP-Ab-02453
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	ABHD14A
Protein Name	Abhydrolase domain-containing protein 14A
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD14A. AA range:201-250
Specificity	ABHD14A Polyclonal Antibody detects endogenous levels of ABHD14A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABHD14A; Abhydrolase domain-containing protein 14A
Observed Band	30kD
Cell Pathway	Cytoplasm . Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Brain,Lymph,
Function	caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,
Background	caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from RAW264.7 cells, using ABHD14A Antibody. The lane on the right is blocked with the synthesized peptide.