



# Cytokeratin 18 Polyclonal Antibody

**Catalog No** YP-Ab-03115

**Isotype** IgG

**Reactivity** Human;Mouse;Rat

**Applications** WB;IHC;IF;ELISA

**Gene Name**

**Protein Name**

**Immunogen**

**Specificity**

**Formulation**

**Source**

**Purification**

**Dilution**

**Concentration** 1 mg/ml

**Purity** ≥90%

**Storage Stability** -20°C/1 year

**Synonyms**

**Observed Band**

**Cell Pathway**

**Tissue Specificity**

**Function**

**Background**

**matters needing attention** Avoid repeated freezing and thawing!

**Usage suggestions** This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

---

---

---

---

---

---

---

---