




ANKRD20 Polyclonal Antibody

Catalog No	YP-Ab-03698
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	ANKRD20A3
Protein Name	Ankyrin repeat domain-containing protein 20A3
Immunogen	The antiserum was produced against synthesized peptide derived from human ANKRD20A3. AA range:551-600
Specificity	ANKRD20 Polyclonal Antibody detects endogenous levels of ANKRD20 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKRD20A3; Ankyrin repeat domain-containing protein 20A3
Observed Band	95kD
Cell Pathway	
Tissue Specificity	Testis,
Function	similarity:Contains 5 ANK repeats.,tissue specificity:Widely expressed; not found in kidney, nor in liver.,
Background	similarity:Contains 5 ANK repeats.,tissue specificity:Widely expressed; not found in kidney, nor in liver.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human brain, using ANKRD20A3 Antibody. The picture on the right is blocked with the synthesized peptide.