



## CATF (Cleaved-Ala271) rabbit pAb

<b>Catalog No</b>	YP-Ab-04366
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB; ELISA
<b>Gene Name</b>	CTSF
<b>Protein Name</b>	CATF (Cleaved-Ala271)
<b>Immunogen</b>	Synthesized peptide derived from human CATF (Cleaved-Ala271)
<b>Specificity</b>	This antibody detects endogenous levels of Human,Mouse,Rat CATF (Cleaved-Ala271, protein was cleaved amino acid sequence between 270-271 )
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	Cathepsin F (CATSF;EC 3.4.22.41)
<b>Observed Band</b>	23 55kD
<b>Cell Pathway</b>	Lysosome.
<b>Tissue Specificity</b>	High expression levels in heart, skeletal muscle, brain, testis and ovary; moderate levels in prostate, placenta, liver and colon; and no detectable expression in peripheral leukocytes and thymus.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

