



GGT5 (heavy chain, Cleaved-Gly387) rabbit pAb

Catalog No	YP-Ab-04383
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	GGT5 GGTLA1
Protein Name	GGT5 (heavy chain, Cleaved-Gly387)
Immunogen	Synthesized peptide derived from human GGT5 (heavy chain, Cleaved-Gly387)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat GGT5 (heavy chain, Cleaved-Gly387, protein was cleaved amino acid sequence between 387-388)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gamma-glutamyltransferase 5 (GGT 5;EC 2.3.2.2;Gamma-glutamyl transpeptidase-related enzyme;GGT-rel;Gamma-glutamyltransferase-like activity 1;Gamma-glutamyltranspeptidase 5;Glutathione hydrolase 5;EC 3.4.19.13;Leukotriene-C4 hydrolase;EC 3.4.19.14) [Cleaved into: Gamma-glutamyltransferase 5 heavy chain; Gamma-glutamyltransferase 5 light chain]
Observed Band	42 65kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images