



## NOTCH3 (Cleaved-Asp1628) rabbit pAb

<b>Catalog No</b>	YP-Ab-04396
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB; ELISA
<b>Gene Name</b>	NOTCH3
<b>Protein Name</b>	NOTCH3 (Cleaved-Asp1628)
<b>Immunogen</b>	Synthesized peptide derived from human NOTCH3 (Cleaved-Asp1628)
<b>Specificity</b>	This antibody detects endogenous levels of Human,Mouse,Rat NOTCH3 (Cleaved-Asp1628, protein was cleaved amino acid sequence between 1628-1629 )
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain]
<b>Observed Band</b>	78 260kD
<b>Cell Pathway</b>	Cell membrane ; Single-pass type I membrane protein.; [Notch 3 intracellular domain]: Nucleus. Following proteolytical processing NICD is translocated to the nucleus.
<b>Tissue Specificity</b>	Ubiquitously expressed in fetal and adult tissues.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images