




VSV-G-Tag Polyclonal Antibody

Catalog No	YP-Ab-04402
Isotype	IgG
Reactivity	Species independent
Applications	WB; ELISA
Gene Name	VSV-G Tag
Protein Name	VSV-G Tag
Immunogen	A synthetic peptide corresponding to aa 501-511 of vesicular stomatitis virus glycoprotein (VSV-G). This vector is a pseudo lentiviral vector.
Specificity	Rabbit polyclonal to VSV-G tag
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-5000, ELISA(peptide)1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Vesicular stomatitis virus glycoprotein tag SV epitope tag,VSV-G epitope tag,vsvg tag
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	Vesicular stomatitis virus (VSV), an enveloped RNA virus from the Rhabdoviridae family, is released from the plasma membrane of host cells by a process called budding. The glycoprotein (VSV-g) contains a domain in its extracellular membrane proximal stem that appears to be needed for efficient VSV budding.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of VSV-G-Tag fusion protein at 0.05-0.15ug, 4° over night, secondary antibody(cat: RS0002) was diluted at 1:10000, 37° 1hour.