



# SARS-E2 Monoclonal Antibody

<b>Catalog No</b>	YP-Ab-04457
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human virus
<b>Applications</b>	WB;ELISA
<b>Gene Name</b>	
<b>Protein Name</b>	
<b>Immunogen</b>	Purified recombinant fragment of SARS-E2 glycoprotein precursor expressed in E. Coli.
<b>Specificity</b>	SARS-E2 Monoclonal Antibody detects SARS-E2 protein.
<b>Formulation</b>	Subclonal supernatant.
<b>Source</b>	Monoclonal, Mouse
<b>Purification</b>	Affinity purification
<b>Dilution</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	
<b>Tissue Specificity</b>	
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western Blot analysis using SARS-E2 Monoclonal Antibody against SARS-E2GP3 recombinant protein.