



SAAL1 Polyclonal Antibody

Catalog No	YP-Ab-05993
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SAAL1
Protein Name	Protein SAAL1
Immunogen	Synthesized peptide derived from human protein . at AA range: 320-400
Specificity	SAAL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	52kD
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in testis and ovary, and to a lesser extent in the lung, spleen and the heart (at protein level).
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the SAAL1 family.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the SAAL1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
