



RFTN2 Polyclonal Antibody

Catalog No	YP-Ab-06011
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	RFTN2 C2orf11
Protein Name	Raftlin-2 (Raft-linking protein 2)
Immunogen	Synthesized peptide derived from human protein . at AA range: 370-450
Specificity	RFTN2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	55kD
Cell Pathway	Cell membrane ; Lipid-anchor .
Tissue Specificity	Brain,
Function	similarity:Belongs to the raftlin family.,
Background	similarity:Belongs to the raftlin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

