



RGF1C Polyclonal Antibody

Catalog No	YP-Ab-06031
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	RASGEF1C
Protein Name	Ras-GEF domain-containing family member 1C
Immunogen	Synthesized peptide derived from human protein . at AA range: 100-180
Specificity	RGF1C Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	51kD
Cell Pathway	intracellular,
Tissue Specificity	Brain,Testis,
Function	function:Guanine nucleotide exchange factor (GEF).,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-GEF domain.,
Background	function:Guanine nucleotide exchange factor (GEF).,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-GEF domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

