



UIMC1 Polyclonal Antibody

Catalog No	YP-Ab-06351
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	UIMC1 RAP80 RXRIP110
Protein Name	BRCA1-A complex subunit RAP80 (Receptor-associated protein 80) (Retinoid X receptor-interacting protein 110) (Ubiquitin interaction motif-containing protein 1)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UIMC1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	79kD
Cell Pathway	Nucleus . Localizes at sites of DNA damage at double-strand breaks (DSBs).
Tissue Specificity	Expressed in testis, ovary, thymus and heart. Expressed in germ cells of the testis.
Function	domain:The Abraxas-interacting region (AIR) mediates the interaction with FAM175A/Abraxas.,domain:The UIM-linker region between the 2 UIM repeats determines the selectivity for 'Lys-63'-linked ubiquitin. The length of the linker is important. The linker reduces the flexibility between the UIM repeats and promotes high-affinity and linkage-selective interactions.,function:Ubiquitin-binding protein that specifically recognizes and binds 'Lys-63'-linked ubiquitin. Plays a central role in the BRCA1-A complex by specifically binding 'Lys-63'-linked ubiquitinated histones H2A and H2AX at DNA lesions sites, leading to target the BRCA1-BARD1 heterodimer to sites of DNA damage at double-strand breaks (DSBs). The BRCA1-A complex also possesses deubiquitinase activity that specifically removes 'Lys-63'-linked ubiquitin on histones H2A and H2AX. Also weakly binds monoubiquitin but with much less aff
Background	ubiquitin interaction motif containing 1(UIMC1) Homo sapiens This gene encodes a nuclear protein that interacts with Brca1 (breast cancer 1) in a complex to recognize and repair DNA lesions. This protein binds ubiquitinated lysine 63 of



histone H2A and H2AX. This protein may also function as a repressor of transcription. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images