



DMRTA Polyclonal Antibody

Catalog No	YP-Ab-06505
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	DMRTA1 DMO
Protein Name	Doublesex- and mab-3-related transcription factor A1
Immunogen	Synthesized peptide derived from human protein . at AA range: 170-250
Specificity	DMRTA Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	55kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.
Function	function:May be involved in sexual development.,similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue specificity:Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.,
Background	function:May be involved in sexual development.,similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue specificity:Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images