



HOMEZ Polyclonal Antibody

Catalog No	YP-Ab-06525
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	HOMEZ KIAA1443
Protein Name	Homeobox and leucine zipper protein Homez (Homeodomain leucine zipper-containing factor)
Immunogen	Synthesized peptide derived from human protein . at AA range: 280-360
Specificity	HOMEZ Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous. Strongly expressed in adult testis and kidney as well as fetal lung and kidney.
Function	function:May function as a transcriptional regulator.,similarity:Contains 2 homeobox DNA-binding domains.,subunit:Homodimer or heterodimer (Potential). Interacts with HOXC8.,tissue specificity:Ubiquitous. Strongly expressed in adult testis and kidney as well as fetal lung and kidney.,
Background	function:May function as a transcriptional regulator.,similarity:Contains 2 homeobox DNA-binding domains.,subunit:Homodimer or heterodimer (Potential). Interacts with HOXC8.,tissue specificity:Ubiquitous. Strongly expressed in adult testis and kidney as well as fetal lung and kidney.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images