



KDM5D Polyclonal Antibody

Catalog No	YP-Ab-06641
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	KDM5D HY HYA JARID1D KIAA0234 SMCY
Protein Name	Lysine-specific demethylase 5D (EC 1.14.11.-) (Histocompatibility Y antigen) (H-Y) (Histone demethylase JARID1D) (Jumonji/ARID domain-containing protein 1D) (Protein SmcY)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	KDM5D Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	169kD
Cell Pathway	Nucleus .
Tissue Specificity	Expression is highly down-regulated in metastatic prostate tumors.
Function	alternative products:Additional isoforms seem to exist,cofactor:Ascorbate.,cofactor:Fe(2+).,domain:The JmjC domain is required for enzymatic activity.,function:Histone demethylase that specifically demethylates 'Lys-4' of histone H3, thereby playing a central role in histone code. Does not demethylate histone H3 'Lys-9', H3 'Lys-27', H3 'Lys-36', H3 'Lys-79' or H4 'Lys-20'. Demethylates trimethylated and dimethylated but not monomethylated H3 'Lys-4'. May play a role in spermatogenesis.,similarity:Belongs to the JARID1 histone demethylase family.,similarity:Contains 1 ARID domain.,similarity:Contains 1 JmjC domain.,similarity:Contains 1 JmjN domain.,similarity:Contains 2 PHD-type zinc fingers.,subunit:Interacts with PCGF6.,
Background	This gene encodes a protein containing zinc finger domains. A short peptide derived from this protein is a minor histocompatibility antigen which can lead to graft rejection of male donor cells in a female recipient. Alternative splicing results



in multiple transcript variants. [provided by RefSeq, Apr 2009],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

