



# PSD4 Polyclonal Antibody

<b>Catalog No</b>	YP-Ab-06796
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse;
<b>Applications</b>	WB;ELISA
<b>Gene Name</b>	PSD4 EFA6B TIC
<b>Protein Name</b>	PH and SEC7 domain-containing protein 4 (Exchange factor for ADP-ribosylation factor guanine nucleotide factor 6 B) (Exchange factor for ARF6 B) (Pleckstrin homology and SEC7 domain-containing protein
<b>Immunogen</b>	Synthesized peptide derived from part region of human protein
<b>Specificity</b>	PSD4 Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	116kD
<b>Cell Pathway</b>	Cell membrane . Cell projection, ruffle membrane . In interphase associated with the plasma membrane, in particular with membrane ruffling regions (PubMed:23603394). Accumulates in dynamic actin-rich membrane ruffles and microvilli-like structures (PubMed:12082148). Recruited to membranes via phosphatidylinositol phosphate-binding (Probable). .
<b>Tissue Specificity</b>	Widely expressed. Highest levels of expression are found in placenta, pancreas, spleen, thymus and peripheral blood.
<b>Function</b>	function:Guanine nucleotide exchange factor for ARF6-like protein. Involved in membrane recycling.,similarity:Contains 1 PH domain.,similarity:Contains 1 SEC7 domain.,subcellular location:Accumulates in dynamic actin-rich membrane ruffles and microvilli-like structures.,tissue specificity:Widely expressed. Highest levels of expression are found in placenta, pancreas, spleen, thymus and peripheral blood.,
<b>Background</b>	function:Guanine nucleotide exchange factor for ARF6-like protein. Involved in membrane recycling.,similarity:Contains 1 PH domain.,similarity:Contains 1 SEC7 domain.,subcellular location:Accumulates in dynamic actin-rich membrane ruffles and microvilli-like structures.,tissue specificity:Widely expressed. Highest levels of expression are found in placenta, pancreas, spleen, thymus and peripheral blood.,



**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**

