



R4RL1 Polyclonal Antibody

Catalog No	YP-Ab-07206
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	RTN4RL1 NGRH2 NGRL2
Protein Name	Reticulon-4 receptor-like 1 (Nogo receptor-like 2) (Nogo-66 receptor homolog 2) (Nogo-66 receptor-related protein 3) (NgR3)
Immunogen	Synthesized peptide derived from human protein . at AA range: 270-350
Specificity	R4RL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	48kD
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor . Membrane raft . Perikaryon . Cell projection . Localized to the surface of neurons, including axons. .
Tissue Specificity	Predominantly expressed in brain. Expressed at lower levels in kidney, lung, mammary gland, placenta, salivary gland, skeletal muscle and spleen.
Function	function:May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.,similarity:Belongs to the Nogo receptor family.,similarity:Contains 8 LRR (leucine-rich) repeats.,subcellular location:Localized to the surface of neurons, including axons.,tissue specificity:Predominantly expressed in brain. Expressed at lower levels in kidney, lung, mammary gland, placenta, salivary gland, skeletal muscle and spleen.,
Background	function:May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.,similarity:Belongs to the Nogo receptor family.,similarity:Contains 8 LRR (leucine-rich) repeats.,subcellular location:Localized to the surface of neurons, including axons.,tissue specificity:Predominantly expressed in brain. Expressed at lower levels in kidney, lung, mammary gland, placenta, salivary gland, skeletal muscle and spleen.,



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images