




F181A rabbit pAb

Catalog No	YP-Ab-07938
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM181A C14orf152 HSD-31 HSD31
Protein Name	F181A
Immunogen	Synthesized peptide derived from human F181A AA range: 42-92
Specificity	This antibody detects endogenous levels of F181A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.52% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein FAM181A
Observed Band	38kD
Cell Pathway	
Tissue Specificity	Hippocampus,Ovary,Testis,
Function	similarity:Belongs to the FAM181 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night