




## DCD2C rabbit pAb

<b>Catalog No</b>	YP-Ab-07951
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human; Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	DCDC2C
<b>Protein Name</b>	DCD2C
<b>Immunogen</b>	Synthesized peptide derived from human DCD2C AA range: 266-316
<b>Specificity</b>	This antibody detects endogenous levels of DCD2C at Human/Mouse
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.65% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	Doublecortin domain-containing protein 2C
<b>Observed Band</b>	38kD
<b>Cell Pathway</b>	Cell projection, cilium, flagellum . Cytoplasm . Detected along the length of the sperm flagellum and in the cytoplasm of the germ cells. .
<b>Tissue Specificity</b>	Expressed in testis and spermatozoa (at protein level).
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western blot analysis of lysates from PC-12 cells,  
primary antibody was diluted at 1:1000, 4° over night