




ZN461 rabbit pAb

Catalog No	YP-Ab-08000
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF461 GIOT1
Protein Name	ZN461
Immunogen	Synthesized peptide derived from human ZN461 AA range: 283-333
Specificity	This antibody detects endogenous levels of ZN461 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.114% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger protein 461 (Gonadotropin-inducible ovary transcription repressor 1) (GIOT-1)
Observed Band	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed, with highest levels in liver, kidney, pancreas, thymus, and small intestine.
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,tissue specificity:Widely expressed, with highest levels in liver, kidney, pancreas, thymus, and small intestine.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night