




1A1L2 rabbit pAb

Catalog No	YP-Ab-08038
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ACCSL
Protein Name	1A1L2
Immunogen	Synthesized peptide derived from human 1A1L2 AA range: 453-503
Specificity	This antibody detects endogenous levels of 1A1L2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.152% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1-aminocyclopropane-1-carboxylate synthase-like protein 2 (ACC synthase-like protein 2)
Observed Band	60kD
Cell Pathway	
Tissue Specificity	
Function	caution:Similar to plant 1-aminocyclopropane-1-carboxylate synthases but lacks a number of residues which are necessary for activity.,similarity:Belongs to the class-I pyridoxal-phosphate-dependent aminotransferase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night