




## RTF1 rabbit pAb

<b>Catalog No</b>	YP-Ab-08418
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human; Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	RTF1 KIAA0252
<b>Protein Name</b>	RTF1
<b>Immunogen</b>	Synthesized peptide derived from human RTF1 AA range: 228-278
<b>Specificity</b>	This antibody detects endogenous levels of RTF1 at Human/Mouse
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1: 500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Nucleus, nucleoplasm .
<b>Tissue Specificity</b>	
<b>Function</b>	function:Plays a role in transcription-coupled histone modification. Plays a role in regulation of transcription. Required for maximal induction of heat-shock genes.,similarity:Contains 1 Plus3 domain.,subunit:Interacts with the RNA polymerase II complex. May interact with the PAF1 complex.,
<b>Background</b>	This locus may represent a gene involved in regulation of transcription elongation and chromatin remodeling, based on studies of similar proteins in other organisms. The encoded protein may bind single-stranded DNA. [provided by RefSeq, Sep 2010],
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western blot analysis of lysates from HpeG2 cells,  
primary antibody was diluted at 1:1000, 4° over night