


CG050 rabbit pAb

Catalog No	YP-Ab-08655
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	C7orf50 FP15621
Protein Name	CG050
Immunogen	Synthesized peptide derived from human CG050 AA range: 100-150
Specificity	This antibody detects endogenous levels of CG050 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4° over night