



CNTD1 rabbit pAb

Catalog No	YP-Ab-09189
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CNTD1 CNTD
Protein Name	Cyclin N-terminal domain-containing protein 1
Immunogen	Synthesized peptide derived from human CNTD1 AA range: 278-328
Specificity	This antibody detects endogenous levels of human CNTD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Cytoplasm . Chromosome . Shuttles between the nucleus and cytoplasm in a stage-specific manner of prophase I cells. Co-localized at crossover sites with PRR19. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

