



YP-Ab-09255

METL5 rabbit pAb

Catalog No	YP-Ab-09255
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	METTL5 DC3 HSPC133
Protein Name	Methyltransferase-like protein 5 (EC 2.1.1.-)
Immunogen	Synthesized peptide derived from human METL5 AA range: 59-109
Specificity	This antibody detects endogenous levels of human METL5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Cell junction, synapse, presynapse . Cell junction, synapse, postsynapse .
Tissue Specificity	Expressed from very early development (8 post-conceptual weeks) and expression persists through adulthood in multiple substructures of the brain, including the cerebellar cortex, hippocampus, and striatum.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

