



SYNG1 rabbit pAb

Catalog No	YP-Ab-09260
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SYNDIG1 C20orf39 TMEM90B
Protein Name	Synapse differentiation-inducing gene protein 1 (SynDIG1) (Dispanin subfamily C member 2) (DSPC2) (Transmembrane protein 90B)
Immunogen	Synthesized peptide derived from human SYNG1 AA range: 195-245
Specificity	This antibody detects endogenous levels of human SYNG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane; Single-pass type II membrane protein. Early endosome membrane ; Single-pass type II membrane protein . Cell junction, synapse, postsynaptic density membrane . Cell junction, synapse . Cell projection, dendrite . Cell projection, dendritic spine . Shuttles between the cell surface and early endosome membrane. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images