



K1468 rabbit pAb

Catalog No	YP-Ab-09269
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KIAA1468
Protein Name	LisH domain and HEAT repeat-containing protein KIAA1468
Immunogen	Synthesized peptide derived from human K1468 AA range: 525-575
Specificity	This antibody detects endogenous levels of human K1468
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Recycling endosome . Golgi apparatus, trans-Golgi network . Translocated to the trans-Golgi network area in an OSBP-dependent manner (PubMed:29514919). Colocalizes with RAB11A in recycling endosomes (By similarity). Found in a complex composed of RELCH, OSBP1 and RAB11A (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images