



## ST65G rabbit pAb

<b>Catalog No</b>	YP-Ab-09287
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	SUPT7L KIAA0764
<b>Protein Name</b>	STAGA complex 65 subunit gamma (Adenocarcinoma antigen ART1) (SPTF-associated factor 65 gamma) (STAF65gamma) (Suppressor of Ty 7-like)
<b>Immunogen</b>	Synthesized peptide derived from human ST65G AA range: 191-241
<b>Specificity</b>	This antibody detects endogenous levels of human ST65G
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Nucleus .
<b>Tissue Specificity</b>	Expressed at high levels in adenocarcinomas and gliomas and low in esophageal cancers and malignant hematological disease. Also expressed at high level in the thymus, low in peripheral blood mononuclear cells, and lowest in the stomach, small intestine, and skeletal muscle.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images