



HS3S1 rabbit pAb

Catalog No	YP-Ab-09318
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HS3ST1 3OST 3OST1
Protein Name	Heparan sulfate glucosamine 3-O-sulfotransferase 1 (EC 2.8.2.23) (Heparan sulfate D-glucosaminyl 3-O-sulfotransferase 1) (3-OST-1) (Heparan sulfate 3-O-sulfotransferase 1) (h3-OST-1)
Immunogen	Synthesized peptide derived from human HS3S1 AA range: 201-251
Specificity	This antibody detects endogenous levels of human HS3S1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Golgi apparatus lumen .
Tissue Specificity	Highly expressed in the brain and kidney and weakly expressed in the heart, lung and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

