



STML2 rabbit pAb

Catalog No	YP-Ab-09320
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STOML2 SLP2 HSPC108
Protein Name	Stomatin-like protein 2 (SLP-2) (EPB72-like protein 2)
Immunogen	Synthesized peptide derived from human STML2 AA range: 51-101
Specificity	This antibody detects endogenous levels of human STML2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Peripheral membrane protein . Mitochondrion . Mitochondrion inner membrane ; Lipid-anchor . Mitochondrion intermembrane space . Membrane raft . Cytoplasm, cytoskeleton . Behaves as an integral membrane protein of the mitochondrion despite the absence of a detectable transmembrane domain (PubMed:21746876). Also associates with the actin cytoskeleton and membrane rafts in activated T-cells (PubMed:18641330, PubMed:10713127). A minor pool is associated with the plasma membrane and is enriched at the immunological synapse in activated T-cells (PubMed:22623988). .
Tissue Specificity	Ubiquitously expressed at low levels. Expressed in lymphoid tissues (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images