



EHD1 rabbit pAb

Catalog No	YP-Ab-09334
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EHD1 PAST PAST1 CDABP0131
Protein Name	EH domain-containing protein 1 (PAST homolog 1) (hPAST1) (Testilin)
Immunogen	Synthesized peptide derived from human EHD1 AA range: 49-99
Specificity	This antibody detects endogenous levels of human EHD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Recycling endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Early endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cell projection, cilium membrane ; Peripheral membrane protein ; Cytoplasmic side . Preferentially associates with tubular recycling endosomes (PubMed:15020713, PubMed:17233914, PubMed:19864458, PubMed:23596323). Colocalizes with FER1L5 at plasma membrane in myoblasts and myotubes (By similarity). Localizes to the ciliary pocket from where the cilium protrudes (PubMed:25686250). Colocalizes with BACE1 in tubulovesicular cytoplasmic membranes. Colocalizes with BACE1 and APP amyloid beta proteins in hippocampal mossy fiber terminals (By similarity). .
Tissue Specificity	Highly expressed in testis.
Function	
Background	



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

