



PK10 rabbit pAbV

Catalog No	YP-Ab-09429
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	
Protein Name	HERV-K_5q33.3 provirus ancestral Pro protein (EC 3.4.23.50) (HERV-K10 Pro protein) (HERV-K107 Pro protein) (Protease) (Proteinase) (PR)
Immunogen	Synthesized peptide derived from human VPK10 AA range: 23-73
Specificity	This antibody detects endogenous levels of human VPK10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
