



WHAMM rabbit pAb

Catalog No	YP-Ab-09555
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	WHAMM KIAA1971 WHDC1
Protein Name	WASP homolog-associated protein with actin, membranes and microtubules (WASP protein homology region 2 domain-containing protein 1) (WH2 domain-containing protein 1)
Immunogen	Synthesized peptide derived from human WHAMM AA range: 387-437
Specificity	This antibody detects endogenous levels of human WHAMM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Endoplasmic reticulum-Golgi intermediate compartment . Cytoplasmic vesicle membrane . Golgi apparatus, cis-Golgi network . Localized to a perinuclear compartment near the microtubule-organizing center (MTOC). Also detected on tubulo-vesicular structures in the cell periphery that frequently localized along microtubules. .
Tissue Specificity	Expressed in brain, lung, heart, colon and kidney (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images