



MOV10 rabbit pAb

Catalog No	YP-Ab-09621
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MOV10 KIAA1631
Protein Name	Putative helicase MOV-10 (EC 3.6.4.13) (Moloney leukemia virus 10 protein)
Immunogen	Synthesized peptide derived from human MOV10 AA range: 13-63
Specificity	This antibody detects endogenous levels of human MOV10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm, P-body . Cytoplasm, Cytoplasmic ribonucleoprotein granule . Cytoplasm, Stress granule . Nucleus . Cytoplasm . Co-enriched in cytoplasmic foci with TUT4 (PubMed:30122351). In developing neurons, localizes both in nucleus and cytoplasm, but in the adulthood it is only cytoplasmic (By similarity). After infection, relocalizes to the DENV replication complex in perinuclear regions (PubMed:27974568). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images