



GNPTG rabbit pAb

Catalog No	YP-Ab-09635
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GNPTG C16orf27 GNPTAG CAB56184 LP2537
Protein Name	N-acetylglucosamine-1-phosphotransferase subunit gamma (GlcNAc-1-phosphotransferase subunit gamma) (UDP-N-acetylglucosamine-1-phosphotransferase subunit gamma)
Immunogen	Synthesized peptide derived from human GNPTG AA range: 216-266
Specificity	This antibody detects endogenous levels of human GNPTG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted . Golgi apparatus .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
