



D107A rabbit pAb

| | |
|----------------------------------|---|
| Catalog No | YP-Ab-09670 |
| Isotype | IgG |
| Reactivity | Human;Rat;Mouse; |
| Applications | WB |
| Gene Name | DEFB107A DEFB107 DEFB7; DEFB107B |
| Protein Name | Beta-defensin 107 (Beta-defensin 7) (BD-7) (DEFB-7) (Defensin, beta 107) |
| Immunogen | Synthesized peptide derived from human D107A AA range: 3-53 |
| Specificity | This antibody detects endogenous levels of human D107A |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution | WB 1:1000-2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | |
| Cell Pathway | Secreted. |
| Tissue Specificity | Specifically expressed in testis. |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images

