



CLPB rabbit pAb

Catalog No	YP-Ab-09704
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CLPB HSP78 SKD3
Protein Name	Caseinolytic peptidase B protein homolog (Suppressor of potassium transport defect 3)
Immunogen	Synthesized peptide derived from human CLPB AA range: 524-574
Specificity	This antibody detects endogenous levels of human CLPB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion intermembrane space .
Tissue Specificity	Widely expressed (at protein level) (PubMed:25597511). Expressed in fetal, as well as in adult tissues, with highest levels in adult brain, including thalamus, hippocampus, occipital cortex and parietal cortex. Low expression in granulocytes (PubMed:25597510).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images