



# TNNI1 rabbit pAb

<b>Catalog No</b>	YP-Ab-09816
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB
<b>Gene Name</b>	TNNI1
<b>Protein Name</b>	Troponin I, slow skeletal muscle (Troponin I, slow-twitch isoform)
<b>Immunogen</b>	Synthesized peptide derived from human TNNI1 AA range: 47-97
<b>Specificity</b>	This antibody detects endogenous levels of human TNNI1
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	
<b>Tissue Specificity</b>	Highest levels observed in human skeletal muscle (e.g. gastrocnemius muscle), differentiated cultures of primary human muscle cells and rhabdomyosarcoma cells cultured in low serum medium. Expressed in C2 muscle cell myoblasts and myotubes.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images