



ME15P rabbit pAb

Catalog No	YP-Ab-09821
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	METTLL15P1 METT5D2
Protein Name	Putative methyltransferase-like protein 15P1 (EC 2.1.1.-) (Methyltransferase 5 domain-containing protein 2) (Methyltransferase-like protein 15 pseudogene 1)
Immunogen	Synthesized peptide derived from human ME15P AA range: 66-116
Specificity	This antibody detects endogenous levels of human ME15P
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

