



CTGE1 rabbit pAb

Catalog No	YP-Ab-09836
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CTAGE1 CTAGE2
Protein Name	Cutaneous T-cell lymphoma-associated antigen 1 (Protein cTAGE-1) (Cancer/testis antigen 21.1) (CT21.1)
Immunogen	Synthesized peptide derived from human CTGE1 AA range: 19-69
Specificity	This antibody detects endogenous levels of human CTGE1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Testis. Expressed in several cutaneous T-cell lymphoma (CTCL) cell lines, in head-neck carcinomas and carcinoma of ovarian tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

