



## SSPN rabbit pAb

<b>Catalog No</b>	YP-Ab-09866
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	SSPN KRAG
<b>Protein Name</b>	Sarcospan (K-ras oncogene-associated protein) (Kirsten-ras-associated protein)
<b>Immunogen</b>	Synthesized peptide derived from human SSPN AA range: 154-204
<b>Specificity</b>	This antibody detects endogenous levels of human SSPN
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cell membrane ; Multi-pass membrane protein . Cell membrane, sarcolemma . Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein . Also found in myotendinous junctions and in the postsynaptic membrane of neuromuscular junctions. .
<b>Tissue Specificity</b>	Isoform 1 is expressed exclusively in heart and skeletal muscle. Isoform 2 is expressed in heart, skeletal muscle, thymus, prostate, testis, ovary, small intestine, colon and spleen.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images