



DOCK9 rabbit pAb

Catalog No	YP-Ab-09894
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DOCK9 KIAA1058
Protein Name	Dedicator of cytokinesis protein 9 (Cdc42 guanine nucleotide exchange factor zizimin-1)
Immunogen	Synthesized peptide derived from human DOCK9 AA range: 1516-1566
Specificity	This antibody detects endogenous levels of human DOCK9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endomembrane system . Associated with membranes. .
Tissue Specificity	Widely expressed, with highest expression in heart and placenta. Expressed at intermediate level in kidney, brain, lung and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
