



VKORL rabbit pAb

Catalog No	YP-Ab-09898
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VKORC1L1
Protein Name	Vitamin K epoxide reductase complex subunit 1-like protein 1 (VKORC1-like protein 1)
Immunogen	Synthesized peptide derived from human VKORL AA range: 28-78
Specificity	This antibody detects endogenous levels of human VKORL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

