



MOS1 rabbit pAb

Catalog No	YP-Ab-09915
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MINOS1 C1orf151
Protein Name	Mitochondrial inner membrane organizing system protein 1
Immunogen	Synthesized peptide derived from human MOS1 AA range: 9-59
Specificity	This antibody detects endogenous levels of human MOS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein . The C-terminus is located in the intermembrane space (By similarity), while the location of the N-terminus has not been determined yet. As some programs predict the presence of 2 closely apposed membrane domains, it has been proposed that the protein may cross the membrane twice and that both termini may face the intermembrane space (PubMed:22114354). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images