



RN19B rabbit pAb

Catalog No	YP-Ab-09943
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RNF19B IBRDC3 NKLAM
Protein Name	E3 ubiquitin-protein ligase RNF19B (EC 6.3.2.-) (IBR domain-containing protein 3) (Natural killer lytic-associated molecule) (RING finger protein 19B)
Immunogen	Synthesized peptide derived from human RN19B AA range: 272-322
Specificity	This antibody detects endogenous levels of human RN19B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic granule membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed specifically in natural killer cells, activated macrophages and cytotoxic T-cells (PubMed:10438909). Present in natural killer cells (at protein level) (PubMed:10912506). Ubiquitously expressed with high expression in testis (PubMed:27485036).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images