



PRCD rabbit pAb

Catalog No	YP-Ab-09982
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PRCD
Protein Name	Progressive rod-cone degeneration protein
Immunogen	Synthesized peptide derived from human PRCD AA range: 4-54
Specificity	This antibody detects endogenous levels of human PRCD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell projection, cilium, photoreceptor outer segment . Membrane ; Lipid-anchor ; Cytoplasmic side . Endoplasmic reticulum . Golgi apparatus . Localizes to photoreceptor disk membranes in the photoreceptor outer segment (PubMed:27613864). The secretion in media described in PubMed:24992209 is probably an experimental artifact (PubMed:24992209). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images