



ALKB8 rabbit pAb

Catalog No	YP-Ab-09754
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ALKBH8 ABH8
Protein Name	Alkylated DNA repair protein alkB homolog 8 (EC 1.14.11.-) (Probable alpha-ketoglutarate-dependent dioxygenase ABH8) (S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8) (tRNA (carboxymethyl
Immunogen	Synthesized peptide derived from human ALKB8 AA range: 27-77
Specificity	This antibody detects endogenous levels of human ALKB8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Nucleus . Predominantly cytoplasmic. .
Tissue Specificity	Widely expressed, with highest expression in spleen, followed by pancreas and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images