



PAR3L rabbit pAb

Catalog No	YP-Ab-09755
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PARD3B ALS2CR19 PAR3B PAR3L
Protein Name	Partitioning defective 3 homolog B (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 19 protein) (PAR3-beta) (Partitioning defective 3-like protein) (PAR3-L protein)
Immunogen	Synthesized peptide derived from human PAR3L AA range: 587-637
Specificity	This antibody detects endogenous levels of human PAR3L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endomembrane system. Cell junction. Cell junction, tight junction. Partially localized along the cell-cell contact region. Colocalizes with TJP1 to epithelial tight junctions.
Tissue Specificity	Highly expressed in kidney, lung and skeletal muscle. Expressed at intermediate levels in brain, heart, placenta, liver and pancreas. Isoform 1 is predominant, while isoform 2 and isoform 3 are expressed at lower levels.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images