



ODR4 rabbit pAb

Catalog No	YP-Ab-09757
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ODR4 C1orf27 TTG1 TTG1A
Protein Name	Protein odr-4 homolog (hODR-4) (LAG1-interacting protein) (Transactivated by transforming growth factor beta protein 1)
Immunogen	Synthesized peptide derived from human ODR4 AA range: 25-75
Specificity	This antibody detects endogenous levels of human ODR4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

